

Instructions for use are to be used in combination with the specific product related information on each product packaging. The gloves are sold as a bundled unit within the shipping carton. In case this bundled unit is dismantled and products are sold separately, the distributor must ensure that the instructions for use are accompanied by each separate unit.

The gloves are classified as Personal Protective Equipment (PPE) Category III according to PPE Regulation (EU) 2016/425 and have been shown to comply with this regulation through the applicable harmonized European standards. These gloves are designed to provide protection against specific chemicals tested, micro-organisms and particulate radioactive contamination (if applicable). The gloves meet the EN/ISO standards shown on each specific packaging. This PPE is single-use only and to be disposed of after contamination.

**Explanation of standards and pictograms**

Permeation levels are based on breakthrough times (tested acc. EN 16523-1:2015+A1:2018) as follows:

ISO 374-1

Type A / B / C



ABCDEFGHIJKLMNOST

Permeation level acc. EN ISO 374-1:2016 +A1:2018

<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>	<b>6</b>
<b>10</b>	<b>30</b>	<b>60</b>	<b>120</b>	<b>240</b>	<b>480</b>

Minimum breakthrough time in minutes

Test chemicals:

A = Methanol / B = Acetone / C = Acetonitrile / D = Dichloromethane / E = Carbon disulphide / F = Toluene / G = Diethylamine / H = Tetrahydrofuran / I = Ethyl acetate / J = n-Heptane / K = Sodium hydroxide 40 % / L = Sulphuric acid 96 % / M = Nitric acid 65 % / N = Acetic acid 99 % / O = Ammonium hydroxide 25 % / P = Hydrogen peroxide 30 % / S = Hydrofluoric acid 40 % / T = Formaldehyde 37 %

Type A = chemical breakthrough time >30 minutes against at least 6 chemicals from the list

Type B = chemical breakthrough time >30 minutes against at least 3 chemicals from the list

Type C = chemical breakthrough time >10 minutes against at least 1 chemical from the list

EN ISO 374-4:2019

The degradation (in %) indicates the change in puncture resistance of the gloves after exposure to the respective challenge chemical.

This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals. The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only and relates only to the chemical tested. It can be different if the chemical is used in a mixture. It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion, and degradation. When used, protective gloves may provide less resistance to dangerous chemicals due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves. Before usage, inspect the gloves for any defects or imperfections.

ISO 374-5: 2016



Virus

Tested for resistance to penetration according to EN ISO 374-2:2019.  
 Tested for resistance to penetration by blood-borne pathogens according to ASTM F1671/F1671M  
 Resistance to bacteria and fungi - pass  
 Resistance to virus - pass  
 The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.

EN 421:2010



Protection against particulate radioactive contamination.  
 These gloves do not protect against mechanical risks.



PPE is for single-use only and must not be reused.




Before usage read instructions for use carefully.



XXXX = Identification number of notified Body responsible for the EU type examination and supervising ongoing conformity

EN ISO 21420:2020 - Protective gloves – General requirements and test methods

**How to don gloves**




Take a glove from the original box and inspect for damages.



Insert hand and pull cuff over wrist until it fits securely.

**How to doff gloves**



Grasp the outside of one glove in wrist area. Peel away to turn it inside out.



Hold removed glove in other hand. Slide fingers under wrist of remaining glove and remove it by rolling it down the hand and fold into the first glove. Discard removed gloves.

**Precautions for use**

Always check the gloves for possible mechanical damage, e.g. holes or tears, before use. Do not use damaged gloves. Glove length is appropriate to the end use where the risk to the wrist area is minimal.

**Storage instructions**

Keep storage area cool, dry and dust free, avoid ventilation and storage close to photocopy equipment. Protect gloves against ultraviolet light sources, sunlight and oxidizing agents. Store in original packaging in a dry and dark place at 10° to 30°C.

**Ingredients / Hazardous ingredients**

Some gloves might contain ingredients which are known to possibly cause skin irritations or allergic reactions with sensitized persons. Check warning information on specific packaging carefully. Formulation available on request.

**Disposal instructions**

Dispose of the gloves in accordance with the valid regulations for this material. Gloves contaminated with chemical substances must be disposed of in accordance with the regulations for the relevant chemicals.