

S semper
med

velvet



**Soft as velvet
and light as a feather**



Hospital



Dental



**Primary
care**



**Long-term
care**



Veterinary

CONTACT US!
sempermed@harpsglobal.com
www.sempermed.com

HARPS Investment Asia Pte. Ltd.
9 Straits View, #08-10A Marina One West Tower
Singapore 018937, Singapore
Tel.: +65 6274 4861 · Fax: +65 6274 6977

HARPS Europe GmbH
Wiedner Guertel 9-13,
1100 Vienna, Austria
Tel.: +43 1 9346348



sempermed® velvet is a lightweight and soft examination and protective glove. It is suitable for a wide range of applications, as indicated by the icons on the cover page, yet its uses are not limited to these. Offering the perfect latex-free alternative for anyone looking for a glove that combines optimal wearing comfort with essential tear resistance. The soft nitrile formulation of the glove makes it wear like a second skin and is especially suitable for people with sensitive skin.

- **Soft as velvet:** the soft material results in high wearing comfort and a comfortable fit.
- **Your daily choice:** suitable for various application fields as per the intended use.
- **Protects the user against chemicals*.**
- **Protects users and patients from the transmission of infections.** Glove material is resistant to blood and other body fluids, and protects against bacteria, fungi and viruses including SARS-CoV-2, Hepatitis B & C and HIV.

Technical data

Type	non-sterile, powder-free examination glove for single use
Color	violet blue
Size / Overall length as per EN 455-2	XS, S, M, L, XL / median 240 mm
Wall thickness measured in single layer in the palm area	min. 0.05 mm typical 0.06 mm
Barrier performance as per EN 455-1	AQL 1.0
Material	Nitrile Butadiene Rubber (NBR)
Glove surface	finger textured
Glove shape	flat with rolled rim, ambidextrous
Force at break as per EN 455-2	median ≥ 6 N
Shelf life in original packaging if stored according to product specification	5 years
Labeling & conformity to standards	EN ISO 21420, EN ISO 374-1 (Type B), EN ISO 374-2, EN 16523-1, EN ISO 374-4, EN ISO 374-5, EN 455-1, EN 455-2, EN 455-3, EN 455-4, ISO 2859-1, ASTM D6319, ASTM F1671/F1671M

Purpose

Single-use medical device class I as set out in MD Regulation (EU) 2017/745



Single-use protective glove – PPE category III (protection against chemical substances for limited time)* as set out in PPE Regulation (EU) 2016/425



Suitable for food contact as set out in Regulation (EC) 1935/2004 on Food Contact Materials



Packaging

available in 100 and 200 pieces per dispenser box
in dispenser box XS, S, M, L 100 / XL 90 pieces
in dispenser box XS, S, M, L 200 / XL 180 pieces
10 dispenser boxes per shipping carton

Article numbers

XS	100pc	3000015965	200pc	3000015946
S	100pc	3000015966	200pc	3000015947
M	100pc	3000015967	200pc	3000015948
L	100pc	3000015968	200pc	3000015949
XL	90pc	3000015969	180pc	3000015950