

New Look & Specifications

Upgraded Examination Glove Packaging

Since September 2023, HARPS Global is the new proud owner of the sempermed brands.

As we embark on a new chapter as HARPS Global, we are excited to unveil the updated portfolio and packaging for our examination gloves, now marketed under the sempermed umbrella brand.

		<p>Medical examination gloves now under the sempermed umbrella brand, replacing sempercare</p> <p>NEW</p>
		<p>Packaging design to incorporate our new branding</p> <p>NEW</p> <p>Our new progressive logo inspired by the infinity symbol, reflects our long-term focus on quality, reliability, & customer service. Semper [sem.per] = latin; always</p>
<p>Semperit Investments Asia Ptd. Ltd.</p>	<p>HARPS Investment Asia Pte. Ltd.</p>	<p>Manufacturer name and address</p> <p>NEW</p>
<p>Semperit Technische Produkte GmbH</p>	<p>EC REP according to MDR (EU) 2017/745 HARPS Europe GmbH</p>	<p>EC Representative, also serving as the Regulatory Importer for EU markets</p> <p>NEW</p>
<p>LOT YYYY-MM XXXXXX YYYY-MM</p>	<p>UDI (xx) xxxxxxxxxxxxxx LOT xxxxxxx (xx) xxxxxx XXXX-XX-XX (xx) xxxxxx XXXX-XX-XX</p>	<p>Implementation of data matrix printing according to MDR UDI on all packaging levels</p> <p>NEW</p>
<p>Barcode icon</p>	<p>Lightbulb icon For detailed product-specific information, please refer to the following page(s).</p>	<p>Article numbers, barcodes, and partial UPDATE of packaging dimensions</p> <p>NEW</p>
<p>5-6 XS 6-7 S 7-8 M 8-9 L 9-10 XL</p>	<p>5-6 XS 6-7 S 7-8 M 8-9 L 9-10 XL</p>	<p>Updated color for size indication</p> <p>UPDATE</p>
<p>AQL ≤ 1.5</p>	<p>AQL AQL 1.0</p>	<p>Change of AQL for pinholes in the core nitrile examination glove portfolio (please refer to individual product specifications for more information)</p> <p>UPDATE</p>
		<p>Enhanced size labeling on outer carton for optimized logistics</p> <p>UPDATE</p>
<p>The product contains dithiocarbamates – not to be used in the case of hypersensitivity to these substances.</p>	<p>Allergen IV</p>	<p>Including the latest symbol from EN 455-3:2023 (replacing the previous accelerator warning)</p> <p>UPDATE</p>
<p>Medical Device Regulation (EU) 2017/745 Personal Protective Equipment (PPE) Cat. III Personal Protective Equipment (PPE) Cat. I </p>	<p>Same conformity to regulations and product classification / same Basic UDI-DIs (please refer to individual product specifications and Declarations of Conformity)</p>	

sempermed® velvet

Great things now come in sempermed packaging!



Our trusted light-weight nitrile examination gloves are now branded as **sempermed® velvet** and will be exclusively available with an AQL 1.0 for pinholes. Please refer to the following reference numbers, GTIN codes:

Reference number	Article Description	GTIN code Dispenser	GTIN-code Carton	Packing	Dispenser Dimensions	Carton Dimensions	Load ability
sempermed® velvet (200) <i>previously sempercure velvet or sempercure nitrile skin² TD</i>							
3000015946	sempermed® velvet 200pc XS	9001570549402	9001570549419				
3000015947	sempermed® velvet 200pc S	9001570549426	9001570549433	200pcs/ dispenser (XS, S, M, L), 180 pcs (XL)	L 240 mm x W 122 mm x H 65 mm	L 340 mm x W 252 mm x H 250 mm	20ft container: 1,392 cartons
3000015948	sempermed® velvet 200pc M	9001570549440	9001570549457				40ft container: 2,838 cartons
3000015949	sempermed® velvet 200pc L	9001570549464	9001570549471	10 dispenser boxes per carton			40ft HC container: 3,174 cartons
3000015950	sempermed® velvet 180pc XL	9001570549488	9001570549495				
sempermed® velvet (100) <i>previously sempercure velvet or sempercure nitrile skin² TD</i>							
3000015965	sempermed® velvet 100pc XS	9001570549754	9001570549747				
3000015966	sempermed® velvet 100pc S	9001570549778	9001570549761	100pcs/ dispenser (XS, S, M, L), 90 pcs (XL)	NEW L 245 mm x W 125 mm x H 42 mm	NEW L 260 mm x W 223 mm x H 255 mm	20ft container: 2,033 cartons
3000015967	sempermed® velvet 100pc M	9001570549785	9001570549792				40ft container: 4,175 cartons
3000015968	sempermed® velvet 100pc L	9001570549815	9001570549808	10 dispenser boxes per carton			40ft HC container: 4,639 cartons
3000015969	sempermed® velvet 90pc XL	9001570549839	9001570549822				

While the **200-piece version** will continue to be offered in the familiar packaging dimensions, the **100-piece version** will feature updated dimensions for improved layering and dispensing.



Hospital



Dental



Primary care



Long-term care



Veterinary



Please contact your Area Sales Manager in case of any questions!

DM3/9_13 Version 1 09/2024 EN

find more information online



CONTACT US!
sempermed@harpsglobal.com
www.sempermed.com

HARPS Investment Asia Pte. Ltd.
9 Straits View, #08-10A Marina One West Tower
Singapore 018937, Singapore
Tel.: +65 6274 4861 · Fax: +65 6274 6977

HARPS Europe GmbH
Wiedner Guertel 9-13,
1100 Vienna, Austria
Tel.: +43 1 9346348