

S semper
med

safe



Strong and durable



Hospital



Dental



Primary
care



Long-term
care



Veterinary

CONTACT US!
sempermed@harpsglobal.com
www.sempermed.com

HARPS Investment Asia Pte. Ltd.
9 Straits View, #08-10A Marina One West Tower
Singapore 018937, Singapore
Tel.: +65 6274 4861 · Fax: +65 6274 6977

HARPS Europe GmbH
Wiedner Guertel 9-13,
1100 Vienna, Austria
Tel.: +43 1 9346348



Feel safe with the strong nitrile glove in Sempermed's examination glove portfolio. Manufactured from nitrile, sempermed® safe offers higher puncture resistance compared to natural rubber gloves with the same thickness and provides longer-lasting chemical resistance against a wide array of industrial chemicals.*

- **The good chemical resistance*** against a wide array of industrial chemicals offers the user enhanced protection and safety.
- Thanks to the absence of natural rubber proteins, the wearer is well-protected against **type I allergies**.
- The textured surface at the glove's fingertips offers a **strong grip**.
- Strong and durable nitrile material makes the glove **suitable for many different uses**.
- **The force at break*** is significantly higher than the standard requirement and therefore provides the user with greater protection.

Technical data

Type	non-sterile, powder-free examination glove for single use
Color	ocean blue
Size / Overall length as per EN 455-2	S, M, L, XL / median 240 mm
Wall thickness measured in single-layer in the palm area	min. 0.08 mm typical 0.09 mm
Barrier performance as per EN 455-1	AQL 1.0
Material	Nitrile Butadiene Rubber (NBR)
Glove surface	finger textured
Glove shape	flat with rolled rim, ambidextrous
Force at break as per EN 455-2	median ≥ 6 N / typical ≥ 9 N
Shelf life in original packaging if stored according to product specification	5 years
Labeling & conformity to standards	EN ISO 21420, EN ISO 374-1 (Type B), EN ISO 374-2, EN 16523-1, EN ISO 374-4, EN ISO 374-5, EN 455-1, EN 455-2, EN 455-3, EN 455-4, EN 421, ISO 2859-1, ASTM D6319, ASTM F1671/F1671M

Purpose

Single-use medical device class I as set out in MD Regulation (EU) 2017/745



Single-use protective glove – PPE category III (protection against chemical substances for limited time)* as set out in PPE Regulation (EU) 2016/425



Suitable for food contact as set out in Regulation (EC) 1935/2004 on Food Contact Materials



Packaging

Dispenser box S, M, L 100 / XL 90 pieces
10 dispenser boxes per shipping carton

Article numbers

S	100pc	3000015996
M	100pc	3000015997
L	100pc	3000015998
XL	90pc	3000015999