

Core Product Portfolio

LEVEL OF PROTECTION (glove weight)

	GOOD	BETTER			BEST	
MAX						
HIGH	 <p>latex edition IC 5.9g XS-XL prev. edition IC</p>	 <p>safe 4.7g S-XL prev. nitrile</p>	 <p>safe+ 6.5g XS-XL prev. safe+</p>			
REGULAR	 <p>skySense +/- 3.2g XS-XL</p>	 <p>latex comfort 5.0g XS-XL prev. edition</p>	 <p>velvet 3.0g XS-XL prev. velvet/nitrile skin²TD</p>	 <p>velvet XP 3.5g XS-XL prev. nitrile skin²</p>	 <p>velvet XP white 3.5g XS-XL</p>	 <p>green X-free 3.3g XS-XL prev. green</p>
LOW	 <p>syn-stretch 5.5g XS-XL prev. syn-stretch</p>	 <p>vinyl comfort 4.8g XS-XL prev. vinyl PF</p>	 <p>latex classic 5.0g XS-XL non EU prev. latex PP</p>			

POSTIONING	LEVEL OF PROTECTION based on glove weight	REF WEIGHT size M	NEW BRAND	CURRENT BRAND	AQL	MATERIAL	SIZES	GLOVE COUNT	MDR	PPE-R	FOOD	GLOVE COLOR
best	regular	3.3g	green X-free	sempercare® green	1.0	nitrile	XS-XL	200 (180 XL)	✓	Cat. III	✓	
better	high	6.5g	safe+	sempercare® safe+	1.0	nitrile	XS-XL	100 (90 XL)	✓	Cat. III	✓	
better	high	4.7g	safe	sempercare® nitrile*	1.0	nitrile	S-XL	100 (90 XL)	✓	Cat. III	✓	
better	regular	3.5g	velvet XP white	n/a	1.0	nitrile	XS-XL	200 (180 XL)	✓	Cat. III	✓	
better	regular	3.5g	velvet XP	sempercare® nitrile skin ²	1.0	nitrile	XS-XL	100 (90 XL) 200 (180 XL)	✓	Cat. III	✓	
better	regular	3.0g	velvet	sempercare® velvet*/nitrile skin ² TD	1.0	nitrile	XS-XL	100 (90 XL) 200 (180 XL)	✓	Cat. III	✓	
good	high	5.9g	latex edition IC	sempercare® edition IC	1.5	latex	XS-XL	100 (90 XL)	✓	Cat. III	✓	
good	regular	3.0-3.5g	skySense	n/a	1.5	nitrile	XS-XL	200 (180 XL)	✓	Cat. III	✓	
good	regular	5.0g	latex comfort	sempercare® edition	1.5	latex	XS-XL	100 (90 XL)	✓	Cat. III	✓	
good	low	5.0g	latex classic (non EU)	sempercare® latex PP	1.5	latex	XS-XL	100 (90 XL)	✓	Cat. I	✓	
good	low	4.8g	vinyl comfort	sempercare® vinyl PF	1.5	vinyl	XS-XL	100 (90 XL)	✓	Cat. I	✓ (for non-fatty food)	
good	low	5.5g	vinyl syn-stretch	sempercare® syn-stretch	1.5	vinyl	XS-XL	100 (90 XL)	✓	Cat. I	✓ (for non-fatty food)	

= Core Portfolio